

Analytical Scale Calibration Data Sheet



Procedure: Document xxxxxx

	Name	Serial #	Calibration	Relevant Information	Manufacturer Repeatability Specs (std dev, measured at, g)
Mass Balance	Mettler Toledo PB3002-S	1126073990	Internal Calibration Fn	Page 11-12 of Manual	Readability: 0.01g, Repeatability (sd): 0.01g, Linearity 0.02, Max Capacity: 3100g
Calibrated Weights	Mettler Toledo	4000022672	Certificate: 01060401-1		

Mettler Toledo PB3002-S

Date: March 27, 2024

Weighing capacity: 25 g/3100 g

Calibratees: Keyur Chuahan, D. Ortiz Rosales

Signature: KC

25 g		
n	Reading (g)	Temp (Deg C)
1	25.01000	
2	25.00000	
3	24.99000	
4	25.01000	
5	24.99000	
6	25.00000	
7	24.99000	
8	25.00000	
9	25.00000	
10	25.01000	
Average	25.0000	
Std. Dev.	0.0082	
% RSD	0.03	
Accuracy (%)	0.00	
Variation (%)	0.0	

10 mg		
n	Reading (g)	Temp (Deg C)
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
Average	#DIV/0!	
Std. Dev.	#DIV/0!	
% RSD	#DIV/0!	
Accuracy (%)	#DIV/0!	
Variation (%)	#DIV/0!	

50 mg		
n	Reading (g)	Temp (Deg C)
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
Average	#DIV/0!	
Std. Dev.	#DIV/0!	
% RSD	#DIV/0!	
Accuracy (%)	#DIV/0!	
Variation (%)	#DIV/0!	

10 mg		
n	Reading (g)	Temp (Deg C)
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
Average	#DIV/0!	
Std. Dev.	#DIV/0!	
% RSD	#DIV/0!	
Accuracy (%)	#DIV/0!	
Variation (%)	#DIV/0!	

501 mg		
n	Reading (g)	Temp (Deg C)
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
Average	#DIV/0!	
Std. Dev.	#DIV/0!	
% RSD	#DIV/0!	
Accuracy (%)	#DIV/0!	
Variation (%)	#DIV/0!	

50 mg		
n	Reading (g)	Temp (Deg C)
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
Average	#DIV/0!	
Std. Dev.	#DIV/0!	
% RSD	#DIV/0!	
Accuracy (%)	#DIV/0!	
Variation (%)	#DIV/0!	

Notes/Observations:

We measure 25g of Uranium oxide at a time to fill the ball mill jars, so we did the initial calibration check using a 20g and 5g calibration weight.