Loc	Target	Pail	Assay 1	A1 date	Assay 2	A2 date	Assay 3	A3 date	Assay 4	A4 date	Assay 5	A5 date	Assayed	Comment
			[mSv/h		[mSv/h		[mSv/h		[mSv/h		[mSv/h]			
1A														
1B														
<b>1C</b>	UC#33	267	39.9	8/6/2021	4.1	1/27/2022							Measured	First one done from the NHC pail changed
2A	UC#36	265	19.6	1/12/2022	10.7	1/28/2022							Measured	0
2B														
2C	C#3	270	37.9	11/10/2021	10.8	1/28/2022							Measured	0
														This pail also contains, wire, a
														wirescanner, IGM guage, target bolts and
3A	ZrC#8	217	205.0	10/2/2019	74.8	12/10/2019	3.1	1/31/2022					Measured	water block bolts
3B	Ta#63	266	410.0	20/01/2022									0	0
3C	Ta#58	211	311.0	8/19/2019	52.5	12/10/2019	7.1	1/31/2022					Measured	0
4A														
4B														
4C														
5A														
5B	SiC#41	260	33.1	13/10/2021	7.8	1/30/2022							Measured	0
5C														
6A	UC#35	264	25.7	11/26/2021	5.9	1/31/2022							Measured	0
6B	Ta#61	245	41.6	1/13/2021	3.6	1/31/2022							Measured	0
6C	Nb#9	246	149.0	1/21/2021	9.1	1/31/2022							Measured	0
														Disc in bottom of pail was not installed pail
														will need repacking. This pail is likely
7A	Ta#62		58.0	9/17/2021		1/27/2022							Measured	contaminated and also contains 2A2M19.
7B	Ta#60	239	83.2	10/21/2020		1/27/2022							Measured	
7C	UC#34	268	419.0	9/22/2021	4.7	1/27/2022							Measured	Pail changed from 258 to 268
8A														
8B														
8C														
_				2111125		211125:-								
9	TM4 opti		1.9	6/14/2018	14.0	3/1/2019	7.1	11/2/2020						In Source Tray Box. The 1.9 mSv/Hr value
	tray 2018													is just for the optics tray. All other values
	TM2 sour													are for both Optics and Source tray.
	tray 2019	)												
10	TM4 Source		22.1	6/10/2019										7.5 to 22.1 mSv/hr on 360 rotation.
	tray													

Need to change pail to new design